# *Biodentine<sup>TM</sup>*

## *Active Biosilicate Technology<sup>TM</sup>*

(FR)	Substitut dentinaire bioactif	1
(EN)	Bioactive dentine substitute	5
(ES)	Sustituto bioactivo de la dentina	9
(DE)	Bioaktiver Dentinersatz	14
(IT)	Sostituto dentinale bioattivo	19
(PT)	Substituto de dentina bioactivo	23
(NL)	Bioactief dentinesubstituut	28
(PL)	Bioaktywny substytut zębiny	32
(EL)	Βιοενεργό υποκατάστατο οδοντίνης	37
(BG)	Биоактивен заместител на дентина	42
(CS)	Biologicky aktivní náhrada dentinu	47
(DA)	Bioaktiv dentinerstatning	51
(ET)	Bioaktiivne dentiini aseaine	55
(FI)	Bioaktiivinen dentiinin korvausaine	59
(HU)	Bioaktív dentinpótló	63
(IS)	Lífvirkt efni til að koma í stað tannbeins	68
(LT)	Biologiškai aktyvus dentino pakaitalas	72
(LV)	Bioloģiski aktīvs dentīna aizstājējs	76
(NO)	Bioaktiv dentinerstatning	80
(RO)	Substitut de dentină bioactiv	84
(SK)	Bioaktívna náhrada dentínu	88
(SL)	Bioaktivni nadomestek dentina	93
(SV)	Bioaktivt dentinsubstitut	97
(HR)	Bioaktivna zamjena za dentin	101

**Utilisation**

Substitut dentinaire bioactif

**Composition**

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Bioactive dentine substitute

#### **Composition**

Powder is composed of tricalcium silicate, zirconium oxide, calcium oxide, calcium carbonate and colourings.

Aqueous solution is composed of calcium chloride and polycarboxylate.

#### **Properties**

Biodentine™ is a bioactive dentine substitute from the « Active Biosilicate Technology™ » innovation.

1. Biodentine™ has mechanical properties similar to the sound dentine and can replace it both in the crown and in the root, without any preliminary conditioning of mineral tissues
2. Biodentine™ contains high-purity, monomer-free mineral ingredients and is highly biocompatible.
3. Biodentine™ creates optimal conditions for maintenance of pulp vitality, by providing a very tight seal on the dentine surface. It therefore reduces **risk of postoperative sensitivity and longevity of restorations in vital teeth.**
4. Biodentine™ being bioactive, creates optimal environment to build reactionary dentine. Dentine bridges are created faster and are thicker than with similar dental materials and represent necessary condition for optimal pulp healing.
5. Biodentine™ offers an initial setting time reduced down to 12 minutes from start of mix, for optimal use in the crown.

#### **Indications**

In the crown:

- Permanent dentine restoration under composites or Inlay/Onlay
- Temporary dentine-enamel restoration.
- Restoration of deep and/or large coronal carious lesions (sandwich technique).
- Restoration of cervical radicular lesions.

On the pulp:

For primary teeth (aged >2-12 years), immature permanent teeth (adolescent aged >12-21 years) and mature permanent teeth (adults aged > 21 years)

- Pulp capping (direct and indirect)

Pulpotomy for diagnosed symptoms of reversible pulpitis and irreversible pulpitis where bleeding is controlled within 5 minutes.

In the root:

- Repair of root perforations.
- Repair of furcation perforations.
- Repair of perforating internal resorptions.
- Repair of external resorption
- Apexification (teeth with open apex).
- Root-end filling in endodontic surgery (retrograde filling).

### **Contra-indications**

Allergy to one of the ingredients.

### **Limitations**

- Restoration of large loss of tooth substance subjected to high stresses.
- Esthetic restoration of anterior teeth.
- Treatment of teeth with irreversible pulpitis where bleeding is not controlled within 5 minutes.

### **Side effects**

Allergic reaction

### **Instructions for use (for each indication)**

#### **Biodentine™ mixing instructions**

1. Open a capsule and place it on the white capsule holder.
2. Detach a single-dose container of liquid.
3. Twist cap to open. Be careful that no drop of liquid falls out of the single dose container.
4. Pour 5 drops from the single-dose container into the capsule.
5. Close the capsule. Place the capsule on a mixing device, such as Technomix, SYG-200, Tac 400 (Lineatac), Silamat, Cap-Mix, Rotomix, Ultramat etc., at a speed of 4000-4200 rotations/min.
6. Mix for 30 seconds.
7. Open the capsule and check the material's consistency.
8. If a thicker consistency is preferred, wait for 30 sec to 1 min before checking again. Do not exceed the working time. If after trituration the material has a granular and not a creamy appearance, add an additional drop of liquid and mix in the amalgamator for an additional 10 seconds. In this case, the setting time may be increased.
9. Collect Biodentine™ with the instrument supplied in the box. Depending on the desired application, you may handle Biodentine™ with an amalgam carrier, a spatula or a Root Canal Messing Gun. Rapidly rinse and clean the instruments to remove any residual material.

### **IMMEDIATE ENAMEL RESTORATION:**

Assess pulp vitality by the usual tests.

1. Isolate the tooth with a rubber dam.
2. Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
3. Adapt a matrix around the tooth if a wall is missing.
4. Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions)
5. Insert Biodentine™ in the cavity, so that the volume of missing dentine is replaced by the same volume of Biodentine™ avoiding to trap air bubbles. Flatten the material without excessive pressure and ensure good adaptation to the cavity walls and margins.
6. Wait until the end of the setting time before performing the permanent enamel restoration. Biodentine™ is compatible with all direct crown restoration techniques and particularly with all types of bonding systems.

## **NON-IMMEDIATE ENAMEL RESTORATION**

Assess pulp vitality by the usual tests.

1. Isolate the tooth with a rubber dam.
2. Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
3. Adapt a matrix around the tooth if a wall is missing.
4. Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions)
5. Insert Biodentine™ in the cavity avoiding trapping air bubbles. Ensure good adaptation of the material to the cavity walls and margins. Do not apply excessive pressure on the material.
6. Model the surface of the restoration
7. Wait until the end of the setting time before removing the matrix
8. To optimize the mechanical properties of the material and facilitate removal of the matrix, a varnish can be applied onto the surface of the restoration.
9. Check occlusion.
10. Within one week to six months after placement of Biodentine™, prepare the cavity according to the criteria recommended for the selected restorative material.

The remaining Biodentine™ material can be considered as sound artificial dentine and permanently left in deep areas of the cavity and in areas adjacent to the pulp chamber. Biodentine™ is compatible with all direct or indirect crown restoration techniques (Inlay/Onlay), and particularly with all types of bonding systems.

## **PULP CAPPING (DIRECT AND INDIRECT):**

Assess pulp vitality by the usual tests.

1. Isolate the tooth with a rubber dam.
2. Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
3. Adapt a matrix around the tooth if a wall is missing.
4. If there is bleeding in the pulp, hemostasis must be achieved before applying Biodentine™.
5. Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
6. Place Biodentine™ directly on the exposed pulp avoiding trapping air bubbles. Ensure good adaptation of the material to the cavity walls and margins. Do not apply excessive pressure on the material.
7. Perform the immediate or non-immediate enamel restoration as indicated above.
8. Patients should be followed according to current recommendations.

## **PULPOTOMY:**

Assess pulp vitality by the usual tests. In case of clinical signs and symptoms of irreversible pulpitis, pulpotomy is recommended when bleeding can be controlled in 5 minutes.

1. Isolate the tooth with a rubber dam.
2. Remove the infected dentine with a round bur and/or a hand excavator.
3. Gain access to the pulp chamber and clean out the pulp.
4. If there is bleeding in the pulp, hemostasis must be achieved before applying Biodentine™. If haemostasis cannot be achieved after 5 minutes, further pulp tissue should be removed (partial or full pulpotomy) step by step until a controlled bleeding.

A full coronal pulpotomy can be carried out to the level of the root canal orifices with bleeding arrested.

5. Adapt a matrix around the tooth if a wall is missing.
6. Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
7. Place Biodentine™ directly in the pulp chamber and ensure good adaptation to the cavity walls and margins
8. Model the surface of the restoration.
9. Wait until the end of the setting time of the material before removing the matrix.
10. To optimize the mechanical properties of the material and facilitate removal of the matrix, a varnish can be applied onto the surface of the restoration.
11. Check occlusion.
12. Within one week to six months after placement of Biodentine™, prepare the cavity according to the criteria recommended for the selected restorative material.
13. Patients should be followed according to current recommendations.

The remaining Biodentine™ material can be considered as sound artificial dentine and permanently left

in deep areas of the cavity and in areas adjacent to the pulp chamber. Biodentine™ is compatible with all direct or indirect crown restoration techniques, and particularly with all types of bonding systems.

#### **REPAIR OF ROOT PERFORATIONS:**

1. Isolate the tooth with a rubber dam.
2. Prepare the root canal alternately using suitable endodontic instruments and a solution of sodium hypochlorite.
3. Dry the canal with paper points and use a chlorhexidine solution or a calcium hydroxide paste for disinfection between visits. Tightly seal the access cavity with a temporary cement to protect the temporary filling.
4. At the next visit (usually after one week), place a rubber dam and remove the temporary crown restoration. Clean the canal alternately using a solution of sodium hypochlorite and suitable endodontic instruments. Dry the canal with paper points.
5. Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
6. Dispense Biodentine™ over the perforation site using a suitable instrument.
7. Condense Biodentine™ with a plugger.
8. Take an X-ray to check that the material is correctly positioned.
9. Remove excess material and place a temporary filling.
10. Complete root canal treatment at the next visit according to current recommendations.

#### **REPAIR OF FURCATION PERFORATIONS:**

1. Isolate the tooth with a rubber dam.
2. Rinse the cavity with a solution of sodium hypochlorite to disinfect the area.
3. If there is bleeding, hemostasis must be achieved before applying Biodentine™.
4. Dry the pulp chamber.
5. Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
6. Dispense Biodentine™ and condense. Perforation repair and crown restoration are performed in a single step.
7. Take an X-ray to check that the material is correctly positioned.
8. Remove excess material.
9. At a subsequent visit, if all clinical signs of a successful treatment are present, the possibility of a permanent restoration can be considered.

#### **REPAIR OF PERFORATING INTERNAL RESORPTIONS:**

1. Isolate the tooth with a rubber dam.
2. Prepare the root canal alternately using suitable endodontic instruments and a solution of sodium hypochlorite.
3. Dry the canal with paper points and use a calcium hydroxide paste for disinfection between visits. Tightly seal the access cavity with a temporary cement to protect the temporary filling.
4. At the next visit (usually after one week), place a rubber dam and remove the temporary crown restoration. Clean the canal alternately using a solution of sodium hypochlorite and suitable endodontic instruments. Dry the canal with paper points.
5. Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
6. Dispense Biodentine™ over the resorptive defect using a suitable instrument.
7. Condense Biodentine™ with a plugger.
8. Take an X-ray to check that the material is correctly positioned.
9. Remove excess material and place a temporary filling.
10. Complete root canal treatment at the next visit according to current recommendations.

#### **APEXIFICATION:**

1. Isolate the tooth with a rubber dam.
2. Prepare the root canal alternately using suitable endodontic instruments and a solution of sodium hypochlorite.
3. Dry the canal with paper points and use a calcium hydroxide paste for disinfection between visits. Tightly seal the access cavity with a temporary cement to protect the temporary filling.
4. At the next visit (usually after one week), place a rubber dam and remove the temporary crown restoration. Clean the canal alternately using a solution of sodium hypochlorite and suitable

endodontic instruments. Dry the canal with paper points.

- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Dispense Biodentine™ in the root canal using a suitable instrument.
- Condense Biodentine™ with a plugger.
- Take an X-ray to check that the material is correctly positioned.
- Remove excess material and place a temporary filling.
- Complete root canal treatment at the next visit according to current recommendations.

#### **ROOT-END FILLING IN ENDODONTIC SURGERY:**

- Gain access to the operative site following the current recommendations in endodontic surgery.
- Using a specific ultrasonic tip, prepare a root-end cavity, 3 to 5 mm deep in the apical portion of the root canal.
- Isolate the area. Achieve hemostasis. Dry the cavity with paper points.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Dispense Biodentine™ in the cavity using a suitable instrument. Condense Biodentine™ with a small plugger.
- Remove excess material and clean the surface of the root.
- Take an X-ray to check that the material is correctly positioned.

#### **Warnings and precautions of use**

- Irritating to eyes
- Ensure that the rubber dam is properly placed so as to completely isolate the operating field.
- Water contamination slows the setting of the material. Prevent exposure to water and fluids during the initial setting stage.
- If adding an additional drop of liquid, wait until the end of the setting time before performing the permanent enamel restoration or removing the matrix, if required.
- Single-use product: this product is designed to be used for one single patient. Reusing it would create a risk of contamination.

#### **Storage**

Do not store above 25°C and protect from moisture.

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