

## Dr. Sascha Herbst


Board Certified Endodontist - Deutsche Gesellschaft für Endodontologie und Traumatologie

Since 2024: Deputy Head of the Department of Conservative Dentistry and Periodontology - University Hospital - Ludwig-Maximilians-University Munich - Munich, Germany

2020 - 2023: Deputy Head of the Department of Oral Diagnostics, Digital Health and Health Services Research - Center for Dental and Craniofacial Sciences - Charité - Universitätsmedizin Berlin - Berlin, Germany

2020 - 2023: Co-Head of the Dental Trauma Board Charité - Universitätsmedizin Berlin - Berlin, Germany

2016 - 2020: Assistant Professor in the Department of Operative and Preventive Dentistry, Center for Dental and Craniofacial Sciences - Charité - Universitätsmedizin - Berlin - Berlin, Germany



# Obturation protocol for general practice

## Summary

### Introduction

Recent advances in material sciences allow for straightforward obturation protocols in general practice.

### Methods

Single-cone obturation using a bioactive calcium silicate-based sealer (BioRoot™ RCS; Septodont; Saint-Maur-des-Fossés, France) was used in two cases, which were followed up for up to 16 months.

### Discussion

Practical benefits include the indication for single-cone obturation, easy post-obturation access to clear the cavity with water, and reduced treatment time and postoperative pain.

### Conclusion

BioRoot™ RCS offers high chances of success for endodontic treatments in general practice.



## Introduction

The primary objective of a root canal treatment is to shape, disinfect and obturate the root canal system to prevent further infection. Besides clinical effectiveness, treatment duration and the pain experience are important to patients as well.

One of the most significant advancements in endodontics in recent years has been the introduction of hydraulic calcium silicate (bioceramic) materials for obturation, like the fully synthetic tricalcium silicate sealer BioRoot™ RCS. Bioceramic sealers are known for their excellent biocompatibility and bioactivity, which means they can stimulate bone mineralisation, facilitate the repair of periapical tissues, and consequently promote the healing process.<sup>(1)</sup> Additionally, a recent systematic review and meta-analysis suggests that bioceramic sealers may be associated with

less post-operative pain and discomfort for patients compared to epoxy resin-based sealers, contributing to a better patient experience.<sup>(2)</sup>

BioRoot™ RCS allows for effective sealing due to tag-like structures at the dentine interface when employing the cold single-cone technique.<sup>(3)</sup> As a result, this material can be used in the single-cone technique, reducing treatment time compared to the warm-vertical compaction technique.<sup>(4)</sup>

Taking the above-mentioned aspects into account, BioRoot™ RCS can be an effective material for use in general practice. The following case reports illustrate the effective use of BioRoot™ RCS in daily practice cases which can be effectively solved by general dentists.

## Case report 1: Second upper molar with straight roots

### Clinical signs and symptoms

A 43-year-old female patient presented to the Department for Oral Diagnostics, Digital Health and Health Services Research at Charité - Universitätsmedizin Berlin with moderate pain while chewing in the area of tooth 27. Additionally, the patient reported a partial fracture of the tooth. Clinically, there was no response to cold testing and no mobility, but tenderness to percussion was observed. The periapical radiograph (*Fig. 1*) revealed an extremely deep caries lesion reaching into the inner part of the dentine, but no signs of a periapical lesion were visible.

### Diagnosis

The clinical and radiographic findings indicated an extremely deep caries lesion with concomitant pulp necrosis and symptomatic apical periodontitis.

### Procedure and treatment

At the first visit, non-selective caries removal was performed, and a pre-endodontic build-up was placed using a bulk-fill composite (SDR, Dentsply-Sirona, Charlotte, USA) (*Fig. 2*). Three root canals were successfully located, and the working length was confirmed



Fig. 01



Fig. 02




Fig. 03



Fig. 04a



Fig. 04b




Fig. 04c




Fig. 05a



Fig. 05b



Fig. 05c




Fig. 06

radiographically (*Fig. 3*). Subsequently, all root canals were enlarged with HyFlex CM (Coltène, X) up to 40.04. During the preparation, pus exudation was observed from the palatal root canal. All root canals were thoroughly irrigated with 3% sodium hypochlorite (NaOCl) and temporarily sealed with a calcium hydroxide intracanal dressing and glass ionomer cement.

After ten days, the patient was free of any symptoms. The final irrigation was done with activated 17% ethylenediaminetetraacetic acid (EDTA) and 3% NaOCl. Before obturation, all canals were rinsed with 0.9% saline solution and dried using paper points (*Fig. 4a-c*). BioRoot™ RCS, which is a bioactive calcium silicate root canal sealer, was manually mixed with five drops of the liquid for a duration of 60 seconds on a sterile glass surface. The sealer was applied into the root canals using a gutta-percha point with gentle, oscillating movements, following the manufacturer's guidelines. Subsequently, the properly fitted gutta-percha points were coated with BioRoot™ RCS, inserted into the root canals, and trimmed at the orifice level using a heat-plugger. The upper part of the melted gutta-percha was meticulously condensed to ensure a suitable

coronal seal. Ultimately, any excess material was flushed out with water and dried using compressed air (*Fig. 5a-c*). The quality of the root canal filling was assessed through a postoperative periapical radiograph (*Fig. 6*).

### Follow up

16 months later, the root canal treatment was so far successful as the patient was still free of any symptoms and the radiograph showed no signs of an apical lesion (*Fig. 7*).




Fig. 07



## Case report 2:

# Lower first molar with straightforward anatomy

### Clinical signs and symptoms

A healthy 30-year-old patient visited the clinic with a previously accessed tooth 46. During the first visit, the patient was pain-free, but she described a history of persistent severe pain before beginning root canal treatment. The preoperative radiograph (Fig. 8) shows the access cavity and an apical lesion at the mesial root.

### Diagnosis

By amalgamating the patient's information and the preoperative radiograph, the initial diagnosis strongly suggested symptomatic apical periodontitis accompanied by an apical lesion.

### Procedure and treatment

The treatment was scheduled for two appointments. At the first visit, the initial access cavity (Fig. 9a, b) was refined and the four root canals were located. Then, all root canals were carefully prepared up

to 40.04 (HyFlex CM) under copious irrigation with 3% NaOCl (Fig. 10a-c). The preparation and working length were confirmed by a periapical radiograph (Fig. 11). Calcium hydroxide was placed into the canals, and the cavity was temporarily sealed with glass ionomer cement.

At the second visit, the patient was still symptom-free and all canals were obturated using BioRoot™ RCS, as previously described in Case 1 (Fig. 12a-c). Finally, the assessment of the root canal filling was conducted using a postoperative periapical radiograph (Fig. 13).

### Follow up

At the follow-up visit 12 months later, the patient was pregnant, so a periapical radiograph was not indicated. Nevertheless, the patient remained symptom-free during the observation period.



Fig. 08



Fig. 09a



Fig. 09b




Fig. 10a



Fig. 10b



Fig. 10c



Fig. 11




Fig. 12a



Fig. 12b




Fig. 12c



Fig. 13

## Discussion

In summary, using a bioceramic root canal sealer offers general dentists a range of advantages, including enhanced biocompatibility and bioactivity compared to epoxy resin-based sealers, which display initial cytotoxic effects.<sup>(5)</sup> Bioactive sealers, on the other hand, have a favourable impact on immune cells and the inflammatory process.<sup>(1,6)</sup> Furthermore, the set material fosters periodontal ligament cell adhesion and aids in the healing process by releasing calcium and silicon ions, along with calcium hydroxide.<sup>(7,8)</sup>

From a practical point of view, the composition of BioRoot™ RCS enables straightforward post-obturation cleaning of the access cavity with water. When combined with the single-cone technique, this material can significantly reduce treatment time.

All the factors mentioned can enhance patient outcomes by reducing both treatment time and postoperative pain.<sup>(2,4)</sup>

## Conclusion

The described single-cone obturation technique using BioRoot™ RCS allows for successful treatment and patient outcomes in general practice.

## References

---

1. Jeanneau C, Giraud T, Laurent P, About I. BioRoot RCS Extracts Modulate the Early Mechanisms of Periodontal Inflammation and Regeneration. *Journal of Endodontics*. 2019 Aug;45(8):1016–23.
2. Seron MA, Nunes GP, Ferrisse TM, Strazzi-Sahyon HB, Victorino FR, Dos Santos PH, et al. Postoperative pain after root canal filling with bioceramic sealers: a systematic review and meta-analysis of randomized clinical trials. *Odontology*. 2023 Oct;111(4):793–812.
3. Atmeh AR, Chong EZ, Richard G, Festy F, Watson TF. Dentin-cement interfacial interaction: calcium silicates and polyalkenoates. *J Dent Res*. 2012 May;91(5):454–9.
4. Andréa Cardoso Pereira, Celso Kenji Nishiyama, Lidiane De Castro Pinto. Single-cone obturation technique: a literature review. *RSBO*. 2013 Dec 13;9(4):442–7.
5. Jung S, Libricht V, Sielker S, Hanisch MR, Schäfer E, Dammaschke T. Evaluation of the biocompatibility of root canal sealers on human periodontal ligament cells ex vivo. *Odontology*. 2019 Jan;107(1):54–63.
6. Castro-Jara S, Antilef B, Osbén C, Alcántara R, Fraga M, Nova-Lamperti E, et al. Bioactivity analysis of calcium silicate-based sealers and repair cements on the phenotype and cytokine secretion profile of CD14+ monocytes: An ex vivo study. *Int Endod J*. 2023 Jan;56(1):80–91.
7. Carlisle EM. Silicon: a possible factor in bone calcification. *Science*. 1970 Jan 16;167(3916):279–80.
8. Han L, Okiji T. Uptake of calcium and silicon released from calcium silicate-based endodontic materials into root canal dentine: Calcium silicate-based endodontic materials. *International Endodontic Journal*. 2011 Dec;44(12):1081–7.